# SEM EDS performance parameters

June 1, 2023

## 1 1. Import the packages

```
[]: import hyperspy.api as hs
import numpy as np
import pandas as pd

from helper_functions import elementlines, nearestlines, theoretical_energy
```

## 3. Select spectra, an specify settings like $i_b$ , $E_0$ , ICR, and PT

```
[ ]: | ##### SEM Apreo March 2023 #####
     path = '../../Masteroppgave/2023-03-08_EDS-Apreo/exports/'
     # All data files are available in the repository at:
     # https://github.com/brynjarmorka/eds-sem-bulk-corrections/tree/main/data
     ##### GaAs #####
     elements = ['Ga', 'As']
     lines_of_interest = ['Ga_Ka', 'As_Ka', 'Ga_La', 'As_La', 'Ga_Kb', 'As_Kb']
     fiori_peaks = ['Ga_La', 'As_La', 'Ga_Ka']
     line_ratio_list = [['As_Ka', 'As_La'], ['Ga_Ka', 'Ga_La'], ['Ga_La', 'As_La']]
     # setting file name, beam current, input count rate and process time
     file_current_ICR_PT = [
         ['GaAs_05kV_25pA.emsa', 25, 880, 6], # 0
         ['GaAs_10kV_25pA.emsa', 25, 1750, 6], # 1
         ['GaAs_15kV_25pA.emsa', 25, 3300, 6], # 2
         ['GaAs_30kV_25pA.emsa', 25, 8000, 6], # 3
         ['GaAs_30kV_50pA.emsa', 50, 16400, 6], # 4
     ]
     ##### GaSb #####
     # Sb_Ka has to high energy
     \# elements = ['Ga', 'Sb']
     # lines_of_interest = ['Ga_Ka', 'Ga_Kb', 'Ga_La', 'Sb_La', 'Sb_Lb1']
     # fiori_peaks = ['Ga_La', 'Sb_La', 'Ga_Ka', 'Sb_Lb1']
```

```
# line ratio list = [['Ga_Ka', 'Ga_La'], ['Sb_La', 'Sb_Lb1'], ['Ga_La', ]
      'Sb_La']]
     # file current ICR PT = [
           ['GaSb 05kV 50pA.emsa', 50, 1080, 6], # 0
           ['GaSb 10kV 50pA.emsa', 50, 2300, 6], # 1
           ['GaSb_15kV_50pA.emsa', 50, 5700, 6], # 2
           ['GaSb_15kV_200pA.emsa', 200, 22000, 6], # 3
     #
           ['GaSb_15kV_400pA.emsa', 400, 42000, 6], # 4
           ['GaSb_30kV_50pA.emsa', 50, 17000, 6], # 5
     #
     #
           ['GaSb_30kV_50pA_noPPUC.emsa', 50, 17000, 6], # 6
     #
           ['GaSb_30kV_50pA_processTime1.emsa', 50, 17000, 1], # 7
     #
           ['GaSb_30kV_50pA_processTime2.emsa', 50, 17000, 2], # 8
           ['GaSb_30kV_50pA_processTime4.emsa', 50, 17000, 4], # 9
           ['GaSb_30kV_400pA_processTime1.emsa', 400, 160000, 1], # 10
     # ]
     # change 'number' to change file
     number = 4
     file = file current ICR PT[number][0]
     beam_current = file_current_ICR_PT[number][1]
     input_count_rate = file_current_ICR_PT[number][2]
     process_time = file_current_ICR_PT[number][3]
[]: # common settings
     make_info_on_all_lines = True
     zero_peak_end_index = 40
     model_background_order = 6
```

3 3. Import the data with HyperSpy, set the elements in the spectrum, and slice off the noise peak

# 4 4. Calculate the Duane-Hunt limit, and slice the spectrum to the limit

```
[]: # # Duane-Hunt method to find the real E_O
     def calculate duane hunt(s=s, buffer start=2, buffer end=0.1,,
      →xaxis_plot_buffer=0.5, dh_plot=False):
         x_max = s.axes_manager[0].high_value # highest x-axis value in keV, used_
      ⇔in Duane-Hunt
         if Vacc > x_max:
            print(f'Vacc={Vacc} > x_max={x_max}, Duane-Hunt not possible')
            return np.nan
         else:
             s dh = s.deepcopy()
             # making the lin fit of the background right before Vacc
             dh start = Vacc-buffer start
             dh_end = Vacc-buffer_end
             s_end = s_dh.isig[dh_start:dh_end] # slice with keV
            m_end = s_end.create_model(auto_background=False)
             m_end.add_polynomial_background(order=1)
            m_end.fit()
             x_s_end = s_dh.isig[dh_start-xaxis_plot_buffer:
      →dh_end+xaxis_plot_buffer].axes_manager[0].axis
             dh_bg_zero_index = np.argmin(np.abs(m_end[-1].function(x_s_end) * s_dh.
      ⇔axes_manager[0].scale))
             dh_limit = x_s_end[dh_bg_zero_index]
             print(f'Duane-Hunt limit: {dh_limit:.3f} keV')
             return dh_limit
```

```
[]: dh_limit = calculate_duane_hunt(dh_plot=False)
```

 $Vacc=30.0 > x_max=20.27$ , Duane-Hunt not possible

```
[]: # Using Duane-Hunt to slice the spectrum
def use_dh_to_slice_spectrum(dh_limit=dh_limit, s=s, plot=False):
    if np.isnan(dh_limit):
        print('No Duane-Hunt limit found, not slicing the spectrum')
    else:
        s = s.isig[:dh_limit]
        print(f'Spectrum sliced at {dh_limit:.2f} keV')
        if plot:
            s.plot(xray_lines=True)
    return s
```

```
[]: s = use_dh_to_slice_spectrum(dh_limit, plot=False)
x = s.axes_manager[0].axis # x-axis in keV, after slicing
```

No Duane-Hunt limit found, not slicing the spectrum

### 5. Make a model of the spectrum, and fit it to the data

```
[]: m = make_model(s=s, model_background_order=model_background_order, plot_m=False)
```

```
line_ratio_to_remove = []
for linepair in line_ratio_list:
    if linepair[0] not in lines_in_model or linepair[1] not in_
lines_in_model:
    # print(f'NB!: {linepair} not in model, ratio not possible')
    line_ratio_to_remove.append(linepair)

# remove the lines from the lists
for line in lines_to_remove:
    lines_of_interest.remove(line)
for linepair in line_ratio_to_remove:
    line_ratio_list.remove(linepair)
```

```
[ ]: remove_lines_not_in_model()
```

#### 6 6. Calibrate the offset and scale

```
[]: lines_of_interest = sort_lines(lines_of_interest)
all_lines = [l.name for l in m if not l.isbackground]
all_lines = sort_lines(all_lines)
```

```
[]: def calibrate_axis(m=m, s=s, rounds=2, xray_lines=lines_of_interest):
    """Calibrating the scale and offset of the energy axis."""
    print('Calibrating energy axis (with many elements it can take multiple_
    →minutes)')

for i in range(rounds):
    print(f'Calibrating offset and scale, round {i+1} of {rounds}')
    m.calibrate_energy_axis(calibrate='offset', xray_lines=xray_lines)
    offset_list.append(s.axes_manager[0].offset)
    m.calibrate_energy_axis(calibrate='scale', xray_lines=xray_lines)
    scale_list.append(s.axes_manager[0].scale)
```

```
print(f'Scale: {scale list[-1]:.6f} eV/px \nOffset: {offset_list[-1]:.6f}_\text{\text{\text{offset}}}
      ⇔keV')
         return scale_list[-1], offset_list[-1]
[]: scale, offset = calibrate_axis(rounds=2, xray_lines=lines_of_interest)
     # using all lines gives basically the same result
     # scale2, offset2 = calibrate axis(rounds=1, xray lines=all lines)
    Calibrating energy axis (with many elements it can take multiple minutes)
    Calibrating offset and scale, round 1 of 2
    Calibrating offset and scale, round 2 of 2
    Scale: 0.010046 eV/px
    Offset: 0.193846 keV
        7. Calibrate the energy resolution
[]: def calibrate_resolution(m=m, s=s, rounds=2, xray_lines=lines_of_interest[:3]):
         Calibrating the energy resolution, i.e. the energy resolution at the Mn Ka_{\sqcup}
         NB! The reference line is the first line in xray_lines, thus the first line ∪
      ⇔should be well defined.
         Using the three strongest lines to first calibrate the width, then estimate \Box
      → the energy resolution,
         which is done with only the first line in lines_of_interest
```

```
NB! The reference line is the first line in xray_lines, thus the first line_□ 
⇒ should be well defined.

Using the three strongest lines to first calibrate the width, then estimate_□ 
⇒ the energy resolution, 
which is done with only the first line in lines_of_interest
"""

for i in range(rounds):
    print(f'Calibrating energy resolution, round {i+1} of {rounds}')
    m.calibrate_energy_axis(calibrate='resolution', xray_lines=xray_lines)
    energy_res_list.append(s.metadata.Acquisition_instrument.SEM.Detector.

⇒EDS.energy_resolution_MnKa)

print(f'Calibrated energy resolution: {energy_res_list[-1]:.3f} eV')
    return energy_res_list[-1]
```

```
Calibrating energy resolution, round 1 of 2
Calibrating energy resolution, round 2 of 2
Calibrated energy resolution: 128.676 eV

[]: def print_calibration_info(scls=scale_list, offs=offset_list, offs=energy_res_list):
    # make pretty print of calibration info
```

[]: energy\_resolution = calibrate\_resolution()

#### []: print\_calibration\_info()

	Current	Original	$\Delta$ original	$\Delta$ last step
Scale [eV/channel]	0.010046	0.010000	0.461 %	0.282 %
Offset [keV]	0.193846	0.200000	-3.134 %	-1.293 %
E-res [eV]	128.676	130.000	-1.019 %	0.031 %

# 8 8. Calibrate the energy and width of the peaks

```
[]: calibrate_lines(rounds=2, xray_lines=lines_of_interest)
```

Calibrating peak positions, round 1 of 2 Calibrating peak positions, round 2 of 2  $\,$ 

# 9 9. Calculate Fiori P/B, peak intensities, FWHMs, and peak deviations

```
[]: def make_lines_info(m=m, all_lines=make_info_on_all_lines, sort_by='Area',

lines_of_interest=lines_of_interest):

lines_info = {}

for i in range(len(m) - 1): # last component is the background

if (all_lines == True) or (m[i].name in lines_of_interest):

lines_info[m[i].name] = {
```

```
'Theoretical E [keV]': theoretical_energy(m[i].name),
                'Calibrated E [keV]': np.round(m[i].centre.value, 4),
                'Area': np.round(m[i].A.value, 1),
                'Fiori P/B': np.round(m[i].A.value / (m[-1].function(m[i].

centre.value) * scale), 1),
                'FWHM [eV]': np.round(m[i].fwhm * 1000, 3),
                'Sigma [keV]': np.round(m[i].sigma.value, 4),
                'Height': np.round(m[i].height * scale, 1),
                'FWHM(Mn Ka)' : np.round(np.sqrt(2.5*(5898.7 - m[i].centre.
 \Rightarrowvalue*1000) + (m[i].fwhm*1000)**2), 3),
    lines_info = pd.DataFrame(lines_info).T
    lines_info['Delta E [eV]'] = np.round((lines_info['Calibrated E [keV]'] -__
 →lines_info['Theoretical E [keV]'] ) * 1000, 2)
    lines_info = lines_info.sort_values(by=sort_by, ascending=False)
    return lines_info
lines_info
```

[]: lines\_info = make\_lines\_info(all\_lines=make\_info\_on\_all\_lines, sort\_by='Area')

[]:		Theoretica	l E [keV]	Calibrated	E [keV]	Area	Fio	ri P/B	\
	Ga_La		1.0980		1.1000	332711.7		585.8	
	Ga_Ka		9.2517		9.2537	301654.4		761.5	
	As_Ka		10.5436		10.5456	183007.1		621.5	
	As_La		1.2819		1.2839	140489.7		236.2	
	Ga_Lb1		1.1249		1.1249	55576.2		97.2	
	Ga_Kb		10.2642		10.2662	38822.9		123.6	
	As_Kb		11.7262		11.7282	26698.9		117.4	
	As_Lb1		1.3174		1.3174	23467.4		39.1	
	Ga_Ll		0.9573		0.9573	18099.5		33.1	
	Ga_Lb3		1.1948		1.1948	15338.0		26.4	
	${\tt Ga\_Ln}$		0.9842		0.9842	8347.7		15.2	
	As_Ll		1.1196		1.1196	6924.7		12.1	
	As_Lb3		1.3860		1.3860	6700.0		11.0	
	$As_Ln$		1.1551		1.1551	2710.0		4.7	
		FWHM [eV]	Sigma [keV	] Height	FWHM(Mn	Ka) Del	ta E	[eV]	
	Ga_La	66.866	0.028	4 46960.0	128	.327		2.0	
	Ga_Ka	162.634	0.069	1 17505.2	134	.396		2.0	
	As_Ka	172.548	0.073	3 10009.8	134	.743		2.0	
	As_La	72.009	0.030	6 18412.9	129	.315		2.0	
	Ga_Lb1	67.993	0.028	9 7714.3	128	.676		0.0	
	Ga_Kb	161.035	0.068	2275.3	122	.529		2.0	
	As_Kb	181.136	0.076	9 1391.1	135	.043		2.0	
	As_Lb1	71.444	0.030	3 3100.0	128	.676		0.0	
	Ga_Ll	64.838	0.027	5 2634.5	128	.676		0.0	
	Ga_Lb3	69.266	0.029	4 2089.9	128	.676		0.0	

```
Ga_Ln
           65.355
                         0.0278
                                   1205.5
                                               128.676
                                                                  0.0
As_Ll
           67.895
                                    962.6
                                                                  0.0
                         0.0288
                                               128.676
As_Lb3
           72.634
                         0.0308
                                    870.6
                                               128.676
                                                                  0.0
As_Ln
           68.546
                         0.0291
                                    373.1
                                               128.676
                                                                   0.0
```

### 10 10. Calculate the relevant peak ratios

```
[]: def peak_ratio(line1, line2, m=m):
         # give the K to L ratio of a line, e.g. 'Ga_Ka' to 'Ga_La'
         try:
             m[line1]
             m[line2]
         except ValueError:
             print('line not in model:', line1, line2)
             return np.nan
         return np.round(m[line1].A.value / m[line2].A.value, 3)
     def calculate_all_line_ratios():
         line ratios = {}
         for line_pair in line_ratio_list:
             pair_name = line_pair[0] + '/' + line_pair[1]
             line_ratios[pair_name] = peak_ratio(line_pair[0], line_pair[1])
         line_ratios = pd.DataFrame(line_ratios, index=['Line ratio']).T
         return line_ratios
```

```
[]: line_ratios = calculate_all_line_ratios()
print(line_ratios)
```

```
Line ratio
As_Ka/As_La 1.303
Ga_Ka/Ga_La 0.907
Ga_La/As_La 2.368
```

## 11 11. Save the results in a DataFrame in a ".csv" file

```
def dead_time(s=s, m=m):
    real_time = s.metadata.Acquisition_instrument.SEM.Detector.EDS.real_time
    live_time = s.metadata.Acquisition_instrument.SEM.Detector.EDS.live_time
    dead_time_percent = (real_time - live_time)/real_time*100
    print(f'Dead_time: {dead_time_percent:.1f}%')
    return round(dead_time_percent, 1)
```

```
[]: dead_time_percent = dead_time(s=s, m=m)
```

Dead time: 43.7%

```
[]: def make_output():
         key_output = pd.DataFrame({
             'Name': [name],
             'Nominal beam energy [kV]': [Vacc],
             'Beam current [pA]': [beam_current],
             'Process time' : [process_time],
             'ICR' : [input_count_rate],
             'Dead time [%]': [dead_time_percent],
             'Live time [s]': [round(s.metadata.Acquisition_instrument.SEM.Detector.
      →EDS.live_time, 1)],
             'Duane-Hunt limit [kV]': [dh_limit],
             'Scale [keV]': [scale],
             'Offset [keV]': [offset],
             'Energy resolution [eV]': [energy_resolution],
         })
         for ratio_pair in line_ratios.T.columns:
             key_output[f'Ratio ({ratio_pair})'] = line_ratios.T[ratio_pair][0]
         for line in fiori_peaks:
             try:
                 key_output[f'Fiori PB ({line})'] = lines_info.loc[line, 'Fiori P/B']
             except KeyError:
                 key_output[f'Fiori PB ({line})'] = np.nan
         return key_output
     key_output.T
                                            0
                               GaAs_30kV_50pA
     Name
                                         30.0
    Nominal beam energy [kV]
```

```
[ ]: key_output = make_output()
```

```
[]:
     Beam current [pA]
                                            50
    Process time
                                             6
     ICR
                                         16400
     Dead time [%]
                                          43.7
    Live time [s]
                                         120.0
     Duane-Hunt limit [kV]
                                           NaN
     Scale [keV]
                                      0.010046
     Offset [keV]
                                      0.193846
    Energy resolution [eV]
                                    128.675916
     Total counts
                                     1930388.0
    Background counts
                                      757476.8
    FWTM/FWHM (Ga Ka)
                                         1.829
    FWTM/FWHM (Ga_La)
                                         1.907
    Ratio (As Ka/As La)
                                         1.303
     Ratio (Ga_Ka/Ga_La)
                                         0.907
```

```
Ratio (Ga_La/As_La)
                                        2.368
    Fiori PB (Ga_La)
                                        585.8
    Fiori PB (As_La)
                                        236.2
    Fiori PB (Ga_Ka)
                                        761.5
[]: def save_output():
        # save the output-df to a csv file
        key_output.T.to_csv(f'results/{name}_output.csv')
         # saves the line info df to another csv file
        save_lines_info = True
        if save_lines_info:
            lines_info.to_csv(f'results/lines_info/{name}_lines_info.csv')
[]: save_output()
[]: # ##### snakkes.
     # ![Image](https://folk.ntnu.no/brynjamm/marathon_dabz.gif)
```